


Disposable Ground Pad & Cable



VERTICAL

• Pediatric / Adult also available

Bipolar

Pediatric / Adult

Monopolar

Pediatric / Adult

HORIZONTAL



• Pediatric / Adult also available

Bipolar

Pediatric / Adult


Monopolar

Pediatric / Adult



Clip size
L×W×H=54×31×11 mm

Detachable Ground Pad Cable




Non-Rem Plug

Rem Plug

Hi-Fi Plug



Clip size
L×W×H=66×36×12 mm



Rem Plug




Non-Rem Plug


- Available for autoclave sterilization
- 3m (10ft) Flexible silicone cable

CE-Ready


Ground Pad & Cable




- GM100 (Box/25pcs)
Ground Pad ,Vertical
Pediatric, monopolar
Non-Sterile packing
Aluminium Pouch



- GB100 (Box/25pcs)
Ground Pad ,Vertical
Pediatric, bipolar
Non-Sterile packing
Aluminium Pouch



- GM900 (Box/25pcs)
Ground Pad, Vertical
Adult, monopolar
Non-Sterile packing
Aluminium Pouch




- GB900 (Box/25pcs)
Ground Pad, Vertical
Adult, bipolar
Non-Sterile packing
Aluminium Pouch




Disposable

- GC2VR (Box/50pcs)
Detachable Ground Pad Cable
PVC Cable length: 3 meter
REM plug
Safety lock system
Non-Sterile packing
PE Bag
Clip size : L×W×H=54×31×11mm


Ground Pad & Cable




- GM110 (Box/25pcs)
Ground Pad ,Vertical
Pediatric, monopolar
Cable length: 2.8 meter
Plug type: Non-REM plug
Non-Sterile packing
Aluminium Pouch



- GB110 (Box/25pcs)
Ground Pad ,Vertical
Pediatric, bipolar
Cable length: 2.8 meter
Plug type: REM plug
Non-Sterile packing
Aluminium Pouch



- GM910 (Box/25pcs)
Ground Pad, Vertical
Adult, monopolar
Cable length: 2.8 meter
Plug type: Non-REM plug
Non-Sterile packing
Aluminium Pouch



- GB910 (Box/25pcs)
Ground Pad, Vertical
Adult, bipolar
Cable length: 2.8 meter
Plug type: REM plug
Non-Sterile packing
Aluminium Pouch

Disposable




- GC2VN (Box/50pcs)
Detachable Ground Pad Cable
PVC Cable length: 3 meter
Non-REM plug
Safety lock system
Non-Sterile packing
PE Bag
Clip size : L×W×H=54×31×11mm




Disposable

- GC2HF6 (Box/50pcs)
Detachable Ground Pad Cable
PVC Cable length: 3 meter
6.3mm Hi-Fi plug
Safety lock system
Non-Sterile packing
PE Bag
Clip size : L×W×H=54×31×11mm


Ground Pad & Cable




- GM100H (Box/25pcs)
Ground Pad ,Horizontal
Pediatric, monopolar
Non-Sterile packing
Aluminium Pouch



- GB100H (Box/25pcs)
Ground Pad ,Horizontal
Pediatric, bipolar
Non-Sterile packing
Aluminium Pouch



- GM900H (Box/25pcs)
Ground Pad, Horizontal
Adult, monopolar
Non-Sterile packing
Aluminium Pouch




- GB900H (Box/25pcs)
Ground Pad, Horizontal
Adult, bipolar
Non-Sterile packing
Aluminium Pouch




Reusable

- GC5VR (Box/50pcs)
Detachable Ground Pad Cable
Silicom Cable length: 3 meter
REM plug
Safety lock system
Non-Sterile packing
PE Bag
Clip size : L×W×H=54×31×11mm
50 times (Alcohol Cleaning)
20 times (Auto clave)


Ground Pad & Cable




- GM110H (Box/25pcs)
Ground Pad ,Horizontal
Pediatric, monopolar
Cable length: 2.8 meter
Plug type: Non-REM plug
Non-Sterile packing
Aluminium Pouch



- GB110H (Box/25pcs)
Ground Pad ,Horizontal
Pediatric, bipolar
Cable length: 2.8 meter
Plug type: REM plug
Non-Sterile packing
Aluminium Pouch



- GM910H (Box/25pcs)
Ground Pad, Horizontal
Adult, monopolar
Cable length: 2.8 meter
Plug type: Non-REM plug
Non-Sterile packing
Aluminium Pouch



- GB910H (Box/25pcs)
Ground Pad, Horizontal
Adult, bipolar
Cable length: 2.8 meter
Plug type: REM plug
Non-Sterile packing
Aluminium Pouch



Reusable

- GC5VN (Box/50pcs)
Detachable Ground Pad Cable
Silicom Cable length: 3 meter
non-REM plug
Safety lock system
Non-Sterile packing
PE Bag
Clip size : L×W×H=54×31×11mm
50 times (Alcohol Cleaning)
20 times (Auto clamp)